

# Blood Pressure Monitor

**BUSINESS LINE:** Medical

**GROUP:** Cardiac Monitoring

MEDICAL DEVICE



## Consumables



**CUFFS\***

MEDICAL DEVICE

### PURPOSE

A Sphygmomanometer or Blood Pressure Monitor, is a medical device used to measure, display and record blood pressure and pulse rate. Consisting of a relatively small monitor linked to an inflatable cuff via a plastic tube.

### GOOD TO KNOW

- Different cuff sizes are available, depending on the size of the patient's arm (adult, child, neonatal...)
- Choose the right cuff according to the patient if needed
- A cuff is always included with the device.