

Centrifuge

BUSINESS LINE: Medical
GROUP: Centrifuge

MEDICAL DEVICE



Accessories



ROTOR*

MEDICAL DEVICE



BUCKETS

MEDICAL DEVICE



ADAPTORS

MEDICAL DEVICE

Consumables



TUBE*

MEDICAL DEVICE

PURPOSE

Laboratory centrifuges are motor driven equipment, used to spin liquid samples at high speeds to separate the liquid samples in to their higher and lower densities. This works using the principle of sedimentation.

GOOD TO KNOW

- The Rotors, Buckets and Adapters are selected for the specifics of the study and the tubes needed.
- Centrifuge Calibration is done yearly by sending a newly calibrated centrifuge to the site and returning the original for re-use.

QSAS

DEVICES	ACCESSORIES
<ul style="list-style-type: none"> Centrifuge Hettich, EBA 200 (220V / 50-60Hz) EU 	<ul style="list-style-type: none"> Rotor included Centrifuge Adapters, Hettich, 1054-A (Inventory) Adapter 13x100mm for centrifuge Hettich, 1058 (Inventory)
<ul style="list-style-type: none"> High Speed Centrifuge Hettich, EBA 200 S (220V - 50/60Hz) 	
<ul style="list-style-type: none"> Cooled Centrifuge Hettich, Mikro 220R / ref 2205 (230V / 50Hz) 	<ul style="list-style-type: none"> Fixed Angle Rotor, Hettich, 1015 (Inventory) Centrifuge Adapters, Hettich, 1054-A (Inventory) Adapter 13x100mm for centrifuge Hettich, 1058 (Inventory)
<ul style="list-style-type: none"> Centrifuge Hettich, Universal 320 / ref 1401 (220V - 50/60Hz) 	<ul style="list-style-type: none"> Swing-out Rotor, Hettich, 1554 (Inventory) Bucket, Hettich, 1559 (Inventory) Adapter for up to 5 tubes for centrifuge, Hettich, 1486 (Inventory) Adapter, Hettich, 1482A (Inventory)
<ul style="list-style-type: none"> Refrigerated Centrifuge Hettich, Universal 320R / ref 1406 (220V - 50/60Hz) 	

QINC

DEVICES	ACCESSORIES
<ul style="list-style-type: none"> Centrifuge Hettich, EBA 200 (110V / 50-60Hz) US 	<ul style="list-style-type: none"> Rotor included Centrifuge Adapters, Hettich, 1054-A (Inventory) Adapter 13x100mm for centrifuge Hettich, 1058 (Inventory)
<ul style="list-style-type: none"> High Speed Centrifuge Hettich, EBA 200 S (100 /127V - 50/60Hz) 	
<ul style="list-style-type: none"> Cooled Centrifuge Hettich, Mikro 220R / ref 2205-01 (100-127V / 60Hz) 	<ul style="list-style-type: none"> Fixed Angle Rotor, Hettich, 1015 (Inventory) Centrifuge Adapters, Hettich, 1054-A (Inventory) Adapter 13x100mm for centrifuge Hettich, 1058 (Inventory)
<ul style="list-style-type: none"> Refrigerated Centrifuge Hettich, Universal 320R / ref 1406-01 (120V - 60Hz) 	<ul style="list-style-type: none"> Swing-out Rotor, Hettich, 1554 (Inventory) Bucket, Hettich, 1559 (Inventory) Adapter for up to 5 tubes for centrifuge, Hettich, 1486 (Inventory) Adapter, Hettich, 1482A (Inventory)

QKK

DEVICES	ACCESSORIES
<ul style="list-style-type: none"> Centrifuge Eppendorf, 5702 (100V - 50/60Hz) 	<ul style="list-style-type: none"> Swing-Out Rotor 85mL x 4 places for 5702R Eppendorf, A-4-38 Round Bucket for rotor A-4-38 Eppendorf, 033719 / 5702 722.006 (Pack of 4) Adapters for 4 to 10mL tubes (for Rotor A-4-38), Eppendorf, 5702 735.000 (Pack of 2)[EU ref] Adapters for 2.6 to 7mL tubes (for A-4-38) Eppendorf, 5702719.005
<ul style="list-style-type: none"> Refrigerated Centrifuge Eppendorf, 5702R (100V - 50/60Hz) [Japan] 	