

# Dataloggers/ Thermometers

BUSINESS LINE: LABORATORY  
GROUP: ENVIRONMENT  
MONITORING



## Accessories



PROBE

### PURPOSE

Used to record temperature over long periods with minimum/maximum, programming & alarm functionalities. Dataloggers records the data over a period and then the data can be extracted into a report (usually in PDF format).

Thermometers used to take the temperature in the present moment in time.  
The devices do not record data.

### GOOD TO KNOW

- Calibration of the device is done by Quipment sending a fully calibrated device to replace the original uncalibrated device every 12 months.
- Quipment offers a wide variety of thermometers to meet specific temperature range needs
- A variety of probes are available for specific uses/needs.
- Quipment offers multiple types of different dataloggers from standard to connected devices.
- Quipment can configure the devices according to the user's needs: Alarms, Measurement intervals, Start delays, etc...