

# Product group: Environment Monitoring

BUSINESS LINE: LABORATORY

## The solution that meets all your needs:



### Data Loggers

Used for temperature recording over a period of time. A data logger as Min/Max, programming and Alarm set up functionalities.

A range of probes are available for specific needs.

We offer a variety of data loggers from standard to connected data loggers.

Data loggers can generate detailed reports to show the recorded values.



### Min/Max Thermometers

Used to measure temperature in the present moment in time. A thermometer is not able to record data.

We offer a variety of probes to meet the specific needs of the user.

To note : some data loggers can function as thermometers.



### Weather Monitor

Weather Monitor allow a complete analysis of the state of the atmosphere and climate, including different variables such as temperature, moisture, wind velocity and atmospheric pressure.