

Incubator

BUSINESS LINE: Laboratory
GROUP: Incubator



Devices for verification of calibration



THERMOMETERS / DATA-LOGGERS

PURPOSE

Laboratory incubators are temperature-controlled equipment used for growing cell or microbiological cultures.

Incubators can maintain optimal temperature and other optional parameters such as humidity or CO2 to meet the needs of the user.

GOOD TO KNOW

The verification of calibration is done every year by Quipment sending a data-logger along with instructions to the site to record the temperature of the device and send an extract to Quipment for verification.

A wide range of incubators exist depending on user needs and voltage capabilities on site. Quipment offers different incubators depending on user needs from standard to more advanced for example the Peltier incubators which uses thermoelectric effect for heating and cooling.