

# Infusion Pumps Systems

BUSINESS LINE: Medical  
GROUP: Infusion

MEDICAL DEVICE



## Accessories



POLE CLAMP

MEDICAL DEVICE



DROP SENSOR

MEDICAL DEVICE



IV POLE

MEDICAL DEVICE

## Consumables



EMPTY BAGS\*

MEDICAL DEVICE



IV SETS / LINES / TUBING\*

MEDICAL DEVICE



CATHETERS\*

MEDICAL DEVICE



SYRINGE AND NEEDLE

MEDICAL DEVICE

### PURPOSE

An Infusion pump offers a solution to administrate a controlled dose to the patient from the connected infusion bag via the IV line.

### GOOD TO KNOW

- The calibration verification is done by send a fully calibrated device to replace the original device.
- It is important to indicate the composition of the consumables.