

Spirometer

BUSINESS LINE: MEDICAL
GROUP: RESPIRATORY



MEDICAL DEVICE

Consumables



ADAPTERS

MEDICAL DEVICE



THERMAL PAPER



TURBINE W/ MOUTHPIECE*

MEDICAL DEVICE



NOSE CLIP

MEDICAL DEVICE



BACTERIAL FILTER

MEDICAL DEVICE

PURPOSE

A spirometer is used to measure the volume and speed of air inspired and expired by the lungs, by measuring the ventilation and movement of air inhaled and exhaled.

GOOD TO KNOW

Accessories and consumables are depending on spirometer.