

Product group: Workstation

BUSINESS LINE: LABORATORY

The solution that meets all your needs:



Laminar Flow Hood

Laminar Flow Hoods move air within a designated space that is uniform in both velocity and direction. When combined with HEPA or ULPA filtration, laminar airflows create exceptionally clean air. Air flows can be horizontal or vertical depending on product sensitivity



PCR Workstation / Clean Bench

PCR Workstations are intended for use in non-hazardous applications where biological or biohazard products are not generated, and user protection is not required. This type of hood is used for DNA/RNA applications. The cabin can be decontaminated using UV light.

If you need more information about Workstations, have a look at our dedicated catalog.